

## PRELIMINARY:



EXPERIMENT	OBSERVATION	INFERENCE
Colour of the salt is noted	white/colorless	Absence of $\text{Cu}^{2+}/\text{Ni}^{2+}/\text{Mn}^{2+}$
Odour is noted	a) vinegar like smell	Presence of acetate ( $\text{CH}_3\text{COO}^-$ )
	b) smell of ammonia	Presence of $\text{NH}_4^+$ ion
	c) No characteristic odour	Absence of $\text{NH}_4^+$ or $\text{CH}_3\text{COO}^-$
Solubility is noted	1) Soluble in water	Presence of $\text{NH}_4^+$ <del>and</del> or $\text{Cl}^-$
	2) Insoluble in water	Absence of $\text{NH}_4^+$ or $\text{Cl}^-$

## TEST FOR $\text{Cl}^-$ , $\text{Br}^-$ , $\text{I}^-$

Take a dry test tube, add a pinch of salt into it. Then add 1-2 ml of conc.  $\text{H}_2\text{SO}_4$  to it and heat.

a) Colourless gas which gives a dense white fumes when a glass rod dipped in  $\text{NH}_4\text{OH}$  is shown near the mouth of test tube

Presence of  $\text{Cl}^-$



### Confirmatory Test

To the salt solution add <sup>1-2 ml</sup> dil.  $\text{HNO}_3$ . Then heat it and then cool it. Now add  $\text{AgNO}_3$  to it

white curdy precipitate which dissolves in  $\text{NH}_4\text{OH}$  soln.

Presence of  $\text{Cl}^-$  is confirmed


## TEST FOR ACID RADICAL (Anion)

### Test for $\text{CO}_3^{2-}$ (Carbonate ion)

Take aml of  
dil. HCl in test  
tube then add  
a pinch of  
salt into it.

Brisk  
effervescence  
with the  
liberation  
of colourless,  
odourless  
gas



Presence of  
 $\text{CO}_3^{2-}$



### Confirmatory:

Pass the above  
gas through  
lime water

lime water  
turns milky







M

Page No.:

Date:



or

Dark

No characteristic  
reaction

Absence of  
 $\text{NO}_3^-$

Add  $\text{NH}_4\text{OH}$  drop by drop to  $\text{FeCl}_3 \Rightarrow$  neutral

Page No:  $\text{FeCl}_3$ .

Date:

The red blood coloured soln is divided into two parts.

To one part add dil  $\text{HCl}$ .



Red colour disappears.

Presence of  $\text{CH}_3\text{COO}^-$  confirmed

To the second part add water and boil it


Reddish brown ppt

Presence of  $\text{CH}_3\text{COO}^-$  confirmed.




Cl<sup>-</sup> - greenish  
I<sup>-</sup> - violet  
P<sup>-</sup> - pale yellow

Page No.:
Date:



b) Reddish brown gas      Presence of Br<sup>-</sup>




### Confirmatory Test:

To the salt soln  
add ~~the~~ dil HNO<sub>3</sub>,  
then heat and  
cool it, Then add  
AgNO<sub>3</sub> to it

Pale yellow ppt

Presence of Br<sup>-</sup> is confirmed






Test for SO<sub>4</sub> SULPHATE ( $\text{SO}_4^{2-}$ )

To the salt  
soln add  
 $\text{BaCl}_2$  soln.

white ppt  
insoluble in  
conc.  $\text{HCl}$ .

Presence of  
 $\text{SO}_4^{2-}$

Presence  
of  $\text{CH}_3\text{COO}^-$   
confirmed



Presence  
of  $\text{CH}_3\text{COO}^-$   
confirmed-

Confirmatory


To the salt  
soln add lead  
acetate soln

white ppt  
soluble in  
excess of  
ammonium  
acetate soln

presence  
of  $\text{SO}_4^{2-}$   
confirmed

+  $3\text{Cl}^-$

$\text{Pb} + 2\text{CH}_3\text{COO}^-$   
when



## TEST FOR BASIC RADICAL


### Test for $\text{O}$ group:

To the dry salt  
add conc.  $\text{NaOH}$   
and heat

Smell of  
ammonia

Presence of  
 $\text{O}$  group  
( $\text{NH}_4^+$ )

confirmatory:




confirmatory:

Pass the above  
gas through  
Nessler's reagent

Brown ppt

Presence of  
 $\text{NH}_4^+$  Confir.







## TEST FOR SULPHATE ACETATE

Take salt  
in a watch  
glass add  
equal quantity  
of oxalic acid,  
then add 1-2  
drop of water  
Rub it using  
your finger  
and smell it

Vinegar like  
smell

Presence  
of acetate  
 $\text{CH}_3\text{COO}^-$



Confirmatory


To the salt  
soln add  
neutral  $\text{FeCl}_3$  soln

Red blood  
colouration

To the  
salt soln.  
add potassium  
iodide soln  
KI


Yellow  
ppt  
which  
dissolves in  
hot water  
and reappears  
as golden  
sparkles  
on cooling

Presence  
of  $Pb^{2+}$   
confirmed



X only for viva

Test for group 2



To the salt  
soln, pass  
add dil HCl  
 $H_2S$  gas

Black ppt or  
Yellow ppt



Presence of  
group 2.

## TEST FOR NITRATE ( $\text{NO}_3^-$ )

To the above mixture add paper balls and heat strongly

Plenty of pale brown gas is produced

Presence of  $\text{NO}_3^-$



(pale brown coloured gas)

### Confirmatory Test (Brown ring test)

Take salt soln in the test tube (saturated)  
Add equal volume of freshly prepared  $\text{FeSO}_4$  soln.  
mix it well, then add conc.  $\text{H}_2\text{SO}_4$  along the sides of the test tube without shaking

Brown ring is formed



Presence of  $\text{NO}_3^-$  is confirmed.



X only for viva

c) Violet coloured gas

Presence of  $I^-$ .




Confirmatory:

To the salt solution  
add 1-2 ml of dil.  $HNO_3$   
and heat it, cool it  
Then add  $AgNO_3$ .

dark  
yellow  
ppt




presence of  
 $I^-$  is  
confirmed



or

No characteristic  
reaction

Absence  
of  $Cl^-$ ,  
 $Br^-$ ,  $I^-$



### Blue ash test

to the salt  
soln add  
dil  $\text{HNO}_3$   
and cobalt  
nitrate soln  
then heat  
it, then dip a  
filter paper into  
the soln, burn  
the filter paper

Blue  
ash is  
formed



Presence of  
 $\text{Al}^{3+}$   
confirmed

the soln into  
two parts.

To one part  
add NaOH soln  
dropwise

white ppt


Presence  
of  $Zn^{2+}$   
confirmed



To the second  
part add Potassium  
ferrocyanide


white  
or bluish  
white ppt

Presence  
of  $Zn^{2+}$   
Confirmed





## Test for 4 group




To the salt  
soln add solid  
 $\text{NH}_4\text{Cl}$  and  
add excess  
of  $\text{NH}_4\text{OH}$   
then add  
 $\text{H}_2\text{S}$

dirty  
white ppt


Presence  
of group 4  
( $\text{Zn}^{2+}$ )

### confirmatory



### Confirmatory :

To the  
salt soln  
dissolve the  
white ppt  
in dil  $\text{HCl}$   
boil off  $\text{H}_2\text{S}$   
then add




Test for group 5:

To the salt  
soln add  
solid  $\text{NH}_4\text{Cl}$   
add excess

white ppt

presence  
of group 5

of  $\text{NH}_4\text{OH}$   
then add  
ammonium  
carbonate soln

Confirmatory


redissolve the  
above ppt in  
acetic acid  
then boil off  $\text{CO}_2$



Then divide into  
3 parts  
To the  
first part add  
Potassium Chromate  
( $K_2CrO_4$ )

Yellow  
ppt


Presence of  
 $Ba^{2+}$   
confirm



To the  
second part  
add ammonium  
sulphate soln

White ppt

Presence of  
 $Sr^{2+}$   
Confirmed



To the third  
part add ammonium  
oxalate soln.


White ppt  
while scratching  
the side  
of the test  
tube

Presence  
of  $Ca^{2+}$   
confirmed



conc HCl is used to convert it to chloride bcz it is more volatile.

Date:



Flame test:

~~et~~  $e^-$  transition is <sup>from lower level to higher level and vice versa</sup> the principle

Take a dry watch glass add pinch of salt into it. Then add

2-3 drops of conc. ~~H<sub>2</sub>SO<sub>4</sub>~~ HCl and make

paste of salt. Take the salt in the glass rod. Show it to the flame.

Note the colour of flame.

1) apple green

$\text{Ba}^{2+}$

confirmed

2) crimson red

$\text{Sr}^{2+}$




confirmed

3) Brick red

$\text{Ca}^{2+}$

confirmed

Group reagent




Test for 1 group:

To the salt soln add dil HCl

white ppt

presence of group 1.




Confirmatory:

To the salt soln add Potassium chromate  $\text{K}_2\text{CrO}_4$

beige yellow ppt

Presence of  $\text{Pb}^{2+}$  confirmed





Test for group 3 $Al^{3+}$ ,  $Fe^{3+}$  X viva

To the salt

soln add excess  
of solid  $NH_4Cl$   
then add  
excess of  $NH_4OH$

gelatinous  
white ppt

Presence of  
group 3  
( $Al^{3+}$ )

Confirmatory:

(Blue lake test)

Dissolve the  
above gelatinous  
ppt in dil.  $HCl$ .  
Add blue litmus  
soln. Then add  
 $NH_4OH$  soln  
in excess  
drop by drop

Blue lake  
formation

Presence of  
 $Al^{3+}$   
confirmed




## Test for group 6 ( $Mg^{2+}$ )

To the salt  
soln add solid  
 $NH_4Cl$  and  
then excess of  
 $NH_4OH$ . Then  
add disodium  
hydrogen phosphate  
 $Na_2HPO_4$

white ppt

Presence  
of group  
6.



### Confirmatory

To the salt  
soln add magneson  
reagent

Blueish  
white  
ppt

Presence  
of  $Mg^{2+}$   
confirmed

Bayer's test - test for saturation/unsaturation.

• Take  $\text{KMnO}_4$  sol.  
make it alkaline

→ To the soln of  
organic compound.  
add alkaline  $\text{KMnO}_4$ .

- Violet color disappears

unsat  
- unsat.

- violet color persists

saturated

1) Test for alcohol. (Aliphatic saturated)


Take the soln of  
organic compound,  
dip in blue litmus  
paper.

litmus turns  
red.

Presence of  
OH group

→ Add a piece of Na metal  
to organic compound.

Presence of  
OH group



Confirmatory test

→ Take organic compound,  
then add 2 ml of ethanol




acid into it  
→ Add 1 or 2 ml of  
conc  $H_2SO_4$  and  
heat strongly.

→ Transfer the contents  
into 15-20 ml of 5%  
aq.  $Na_2CO_3$  soln. taken  
in a beaker. Smell the  
content.

pleasant fruity  
smell.

presence  
of OH  
- grp.  
confirmed



2) Test for phenolic grp.

→ Blue litmus test

turns red.

presence  
of phenolic  
grp Ar-OH.


1. To the soln. of organic  
compound add blue litmus

To the aq. soln of organic  
compound add neutral  
ferric chloride soln.

violet coloration.

Presence  
of phenol  
grp.






3) Test for carboxylic acid.  $>\text{C}=\text{O}$   
(Litmus test)

→  $\text{NaHCO}_3$  test

i) To the organic compound  
solid sodium bicarbonate.

Brisk effe-  
-escence with  
the liberation of  
colorless, odourless  
gas



Confirmatory test.

ii) Ester Test

i) To the organic compound.  
with ethanol.  $\text{C}_2\text{H}_5\text{OH}$

4) Test for carbonyl grp.  
2,4 DNP test

→ To the given  
organic compound.  
add 1 ml alcohol

# Organic compounds

papergrid

Date: / /

1. Preliminary Test  
Colour noted Colourless absence of phenol, cresol etc
2. Odour  
Carbolic smell Presence of phenolic grp.  
Pleasant fruity smell Presence of Ester

No characteristic odour Absence of Phenol/Ester

3. Take the given organic compound & burn it Sooty flame aromatic  
non-sooty flame aliphatic

4. Test for unsaturation

Baeyer's test

To the org. compound soln. Violet colour <sup>disappears</sup> unsaturated  
add alk.  $KMnO_4$  Violet colour persists saturated


Test for Alcohol

Take the soln. of org. compound & dip a blue litmus paper. turn red presence of alc. group

blue litmus paper.

5. Reaction with Na

Add a piece of Na metal Evolution of gas Presence of alcoholic group





### Test for ketone

Lower Sodium Nitroprussiate test

Take 1 ml of Sodium nitroprusside soln. and make it alkaline by adding a few drops of NaOH soln.  $\rightarrow$  To it add a pinch of Organic compound and shake it

A red blood coloration

Presence of ketone

### Meta dinitrobenzene test

To the given organic compound add 1 ml of ethanol and then add 1 ml of ethanolic soln of Metadinitrobenzene. and a few drops of NaOH soln. and then shake.


Violet coloration

presence of  $\text{C}=\text{O}$



then add 4 ml of  
alcohol

2, 4 DNP reagent a yellow, orange  
and shake the or red ppt  
mixture. Will separate out



Test for Aldehyde

To the organic compound  
add Schiff's base.

pink color

presence of  
aldehyde

Confirmatory

- i) Tollen's reagent is  
added to organic compound.  
Heat strongly in a water bath.

appearance of  
~~the~~ silver mirror  
or grey ppt

Presence  
of CHO  
grp.

Take equal volume of  
Fehling's reagent no. 1 & 2.  
Then add organic compound  
and heat strongly in  
water bath.

Reddish brown  
ppt.

Prese  
nce of  
-CHO.

### Test for ketone

Lower Sodium Nitroprussiate test

Take 1 ml of Sodium nitroprusside soln. and make it alkaline by adding a few drops of NaOH soln.  $\rightarrow$  To it add a pinch of Organic compound and shake it

A red blood coloration

Presence of ketone

### Meta dinitro benzene test

To the given organic compound add 1 ml of ethanol and then add 1 ml of ethanolic soln of Metadinitro benzene. and a few drops of NaOH soln. and then shake.

Violet coloration

presence of  $\text{t-c-i}$